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
en

Hologic and MyoSure are registered trademarks of Hologic, Inc. and its subsidiaries in the United States and other countries. Aquilex is a trademark of Hologic, Inc. and its subsidiaries in the United States and other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners.

These instructions for use contain information that is subject to copyright. All rights reserved. These instructions for use may not be photocopied, duplicated on microfilm, or otherwise copied or distributed, completely or in part, without the approval of W.O.M. WORLD OF MEDICINE GmbH.

We reserve the right to technical changes without prior notification due to the continuous further development of our products. Function or design may partially differ from the description in the instructions for use. Please contact us for additional information about this or any of our other products.

Some of the parts and equipment referred to in these instructions for use are associated with registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark symbol indicates that any given designation is not subject to trademark protection.

Users of this product should not hesitate to point out to us any errors or issues concerning these instructions for use.

Copyright © W.O.M. WORLD OF MEDICINE GmbH

Federal Law (only for U.S. market)

Caution: Federal law restricts this device to sale by or on the order of a physician.

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4 Initial System Set-up ... 13
4.1 Preparing the System for Use ... 13
4.2 System Components ... 14

5 System Operation ... 15
5.1 Front of Irrigation Pump Unit ... 15
5.2 Rear of Irrigation Pump Unit ... 15
5.3 Fluid Monitoring Unit Set-up ... 16
5.3.1 Setting of the Container Scale ... 18
5.3.2 Connecting the Vacuum Tube ... 19
5.4 Turning On the Aquilex System ... 20
5.5 Hanging the Fluid Bags ... 21
5.6 Using Tube Sets ... 21
5.7 Tube Overview ... 22
5.8 Connecting the Outflow Tube Set ... 23
5.8.1 Connecting Outflow Tube of Tissue Removal Handpiece (e.g. MyoSure®) ... 24
5.9 Inserting the Inflow Tube Set ... 25
5.10 Presetting the Intrauterine Pressure ... 26
5.11 Deficit Limit Setting ... 26
5.12 Using the Pump during Surgery ... 27
5.13 Changing Bags during Surgery ... 28
5.14 Changing Container during Surgery ... 28
5.15 Changing Instrument during Surgery ... 29
5.16 Total Inflow Volume Displayed ... 29
5.17 Turning System Off ... 29

6 Functional Check ... 30

7 Safety Functions ... 31

8 Care and Maintenance ... 32
8.1 Cleaning the System ... 32
8.2 Maintenance Carried out by Authorized Service Technician ... 32
8.3 Replacing of the Fuse ... 33

9 Annual Inspection ... 34
9.1 Safety Test ... 34
9.2 Basic Function Tests ... 34
9.2.1 Scale Test ... 35
9.2.2 Flow Rate Test ... 35
9.2.3 Pressure Measuring Test ... 36
9.2.4 Fluid Deficit Measurement Test ... 37
9.2.5 Testing the Vacuum Pump ... 38
9.3 Determine the Software Version ... 39

10 Error and Warning Messages ... 40

11 Technical Data ... 42

12 Guidelines and Manufacturer's Statement - Electromagnetic Compatibility ... 44
12.1 Electrical Connections ... 44
12.2 Guidelines and Manufacturer’s Statement – Electromagnetic Emissions ... 45
12.3 Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity ... 46
12.4 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity ... 47
12.5 Recommended Safety Distances Between Portable and Mobile RF Telecommunications Devices and Aquilex Fluid Control System ... 48

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1 Important Operator/User Notes

en

Read the manual carefully and become familiar with the operation and function of the Aquilex™ Fluid Control System (Aquilex System) and the accessories before using the device in the OR. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing staff or service personnel, or
- damage or malfunction of the system and/or accessories.

The manual is only for the Aquilex fluid control system, consisting of pump, container scale and bag scale.

Subject to technical changes

The manufacturer reserves the right to have illustrations and technical data deviating slightly from the delivered product due to further development of the products.

Please note

The words WARNING, CAUTION, and NOTE carry special meanings. Read these sections with special attention.

WARNING!
Warnings indicate risks to the safety of the patient or operator. Failure to follow warnings may result in injury to the patient or operator.

CAUTION!
Warnings indicate risks to the equipment. Failure to follow cautions may result in damage to the system.

NOTE!
Notes provide special information to clarify instructions or present additional information.

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2 Safety Information

Exclusion of liability

Hologic is not liable for direct or consequential damage and the warranty is null and void if:

- the system and/or the accessories are operated and used by untrained personnel,
- the system and/or the accessories are improperly used, prepared, or maintained,
- the instructions and rules in the user/operator manual are not adhered to,
- unauthorized persons perform repairs, adjustments, or alterations on or to the system or accessories,
- unauthorized persons open the system,
- the prescribed inspection and maintenance schedule is not adhered to.

Receipt of technical documentation from Hologic does not authorize individuals to perform repairs, adjustments, or alterations on or to the system or accessories.

Authorized service technician

Only an authorized service technician may perform repairs, adjustments, or alterations on the system or accessories. Any violation will void the manufacturer's warranty. Authorized service technicians are trained and certified only by the manufacturer.

Normal Use

The system may be used only as intended.

Care and maintenance

The service and maintenance of the device and its accessories has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is complete and functional before each use.

Maintenance of the device may not be performed during the operation.

Contamination

Before shipping, decontaminate device and accessories in order to protect the service personnel. Follow the instructions listed in this manual. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to reject contaminated products for repair.

Waste management

In the European Community, this symbol indicates that the waste of electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately instead. Please contact Hologic or an accordingly authorized disposal or waste management company for further information.

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3 Purpose

Intended Use

The Aquilex™ Fluid Control System is intended to provide fluid distension of the uterus during diagnostic and operative hysteroscopies and to monitor the volume differential between the irrigation fluid flowing into and out of the uterus.

Contraindications

The system may not be used to introduce fluids into the uterus when hysteroscopy is contraindicated. See the operators manual of your hysteroscope for absolute and relative contraindications.

Relative contraindications to endometrial ablation:

Hysteroscopic endometrial ablation, whether by laser or electrosurgery, should not be undertaken before adequate training, preceptorship, and clinical experience. Additionally, tissue sampling is required prior to destruction of the endometrium. The following are clinical conditions that can significantly complicate hysteroscopic endometrial ablation:

- Adenomatous endometrial hyperplasia
- Uterine leiomyoma
- Severe adenomyosis
- Pelvic pain (subtle PID)
- Uterine anomalies
- Surgical skill (see above)
- Severe anemia
- Inability to circumnavigate the myoma (re: myoma size) - predominantly intramural myomas with small submucous components.

3.1 Warnings and Precautions

3.1.1 Warnings

WARNING!
When performing monopolar hysteroscopic electrosurgery, the distension medium must be electrically non-conductive. Examples include glycine, sorbitol and mannitol. Isotonic saline irrigation fluids may only be used when performing bipolar electrosurgical resective procedures.

WARNING!
The pressure should be kept as low as possible to allow for a sufficient intrauterine distension and to reduce the forces that could allow fluid, ambient air, and/or gas to enter the circulatory system.

WARNING!
Intrauterine distention is usually possible with pressure values between 35 to 70 mmHg. A pressure above 75 to 80 mmHg is required only in rare cases or if the patient has an excessively high blood pressure.

WARNING!
Fluid overload
There is a risk of irrigation fluid reaching the circulatory system of the patient's soft tissue by passing through the uterus. This can be affected by distention pressure, flow rate, perforation of the uterine cavity and duration of the hysteroscopic surgery. It is critical to closely monitor the input and outflow of the distending liquid at all times.

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WARNING!
Fluid deficit
The fluid left in the patient must be monitored. The deficit is the total amount of fluid left in the patient or unaccounted for otherwise. Take notice of the measurement tolerance of the system (see Chapter 11, Technical Data). Estimating the fluid volume remaining in the patient is the physician's responsibility.

WARNING!
Fluid intake and output surveillance
Strict fluid intake and output surveillance should be maintained. If a low viscosity fluid distention medium is used, intrauterine instillation exceeding 2 liters should be followed with great care due to the possibility of fluid overload.

WARNING!
Serum sodium concentration
It is also necessary to monitor the concentration of sodium in the blood of the patient to prevent electrolyte disturbances. Monitoring of the concentration of sodium in the blood must be performed by the physician and is not performed or supported by the system.

WARNING!
The deficit display value is lost in case of a power loss or "brownout."

WARNING!
If the message "Check Scale Connection" appears, the deficit must be calculated manually. The pump keeps displaying the last known deficit value determined prior to the failure of the scale connection.

WARNING!
A container change during surgery is only allowed, if the container holds at least 0,5 liters of fluid. Otherwise, the deficit value may be falsified. In this case, the manufacturer recommends manual deficit calculation.

WARNING!
Hyponatremia
Some distension fluids may lead to fluid overload and, consequently, hyponatremia with its attending sequelae. This can be affected by the distending pressure, flow rate, and duration of hysteroscopic procedure. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!
Pulmonary edema
Hysteroscopic surgery is associated with a risk of developing pulmonary edema resulting from fluid overload with isotonic fluids. It is critical to closely monitor the input and outflow of the distending liquid at all times.

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WARNING!
Cerebral edema
Hysteroscopic surgery is associated with a risk of developing cerebral edema resulting from fluid overload and electrolyte disturbances with hypoosmolar (non-ionic) fluids such as glycine 1.5 % and sorbitol 3.0 %. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!
Idiosyncratic reactions
In rare cases, idiosyncratic reactions, including:
- **intravascular coagulopathy**
- **allergic reaction including anaphylaxis**

may occur while performing a hysteroscopy if a liquid distention medium is used.

WARNING!
Hypothermia (monitoring body temperature)
Continuous flow of distention fluids can lead to a lowering of the patient's body temperature during hysteroscopic surgery. Lower body temperatures can cause coronary and cardiovascular problems. Always monitor the patient's body temperature during the entire surgery. Make especially sure that the following, hypothermia promoting, operation conditions are avoided as best as possible:
- **longer operating times**
- **use of cold irrigation fluid.**

WARNING!
Rupture of the fallopian tube secondary to tubal obstruction
Distention of the uterus may lead to a tear of the fallopian tube should there be an obstruction or permanent occlusion. The rupture could lead to irrigation fluid flowing into the patient's peritoneal cavity, resulting in a fluid overload. It is critical to closely monitor the input and outflow of the distending liquid at all times.

WARNING!
An air embolism can be the result of air contained in the tube set or connected instrument reaching the patient. Ensure there is always fluid in the bag to prevent air from being pumped into the patient.

WARNING!
The system is only intended for use with flexible fluid containers. Do not use glass containers as they might break. With rigid containers, fluid cannot flow quickly enough due to the vacuum being generated inside of the containers. Risk of implosion with rigid containers.

WARNING!
Filling the tubing with irrigation fluid and resetting the deficit display to zero are to be done at the physician’s discretion.

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WARNING!
Place the system in such a way as to allow for easy visualization of the display values, system functions, and access to the control elements.

WARNING!
Functional test
The functional test must be performed prior to each device use.

WARNING!
Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

WARNING!
Pressing the ON/OFF switch does not disconnect the system from the wall power outlet. This requires pulling the power cord located in the rear of the system.

WARNING!
Technique and procedures
Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this system is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical outcome.

WARNING!
Check all factory settings
Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.

WARNING!
Original accessories
For your own safety and that of your patient, use only Aquilex accessories.

WARNING!
Additional equipment
Additional equipment connected to medical electrical devices must be demonstrated to be compliant with their respective IEC or ISO standards (IEC 60601-1, IEC 60950 or IEC 62368 for data processing equipment).

WARNING!
Not explosion-proof
The system is not explosion-proof. Do not use in an area where flammable anesthetic gases are present.

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WARNING!
Changes to the system are not allowed.

WARNING!
To avoid risk of electrical shock, this system may only be connected to a supply mains with protective earth.

WARNING!
Professional qualification
This manual does not include descriptions or instructions for surgical procedures/techniques. It is also not suitable for training physicians in the use of surgical techniques. Medical instruments and systems may be used only by physicians or medical assistants with the appropriate technical/medical qualification working under the direction and supervision of a physician.

WARNING!
Sterile media and accessories
Always work exclusively with sterile substances and media, sterile fluids, and sterile accessories, if so indicated.

WARNING!
Replacement system and accessories
In case the system or any of the accessories fail during a procedure, an alternative system and replacement accessories should be kept within easy reach to be able to finish the operation with the replacement components.

WARNING!
Cleaning the system / Sterilization not allowed.
The pump and the cart/scale can be disinfected by wiping off the outer surfaces. Do not sterilize the pump and the cart/scale.

WARNING!
Condensation / Water penetration
Protect system from moisture. Do not use if moisture has penetrated the system.

WARNING!
System defect
If a system defect is suspected or confirmed, do not use the system. Ensure the system will no longer be used until a qualified service technician conducts the appropriate tests and repairs.

WARNING!
Risk of electrical shock
To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.

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WARNING!
Replacing fuse
Replace the fuse only with a fuse of the same type and rating (see Chapter 11, Technical Data).

WARNING!
Equipment should be positioned such that power cord can be easily disconnected.

WARNING!
Electromagnetic emissions may increase and rise above the permissible limits if other equipment (e.g. MyoSure® Control Unit) is stacked onto or placed directly next to the Aquilex Fluid Control System. The user is responsible for monitoring the devices to make sure they function properly.

WARNING!
Electrical Interference (see Chapter 12, Guidelines and Manufacturer's Statement - Electromagnetic Compatibility): Electrical interference with other devices or instruments was considered when developing this system and none was detected during testing. However, if you still detect or suspect such interference, please follow these suggestions:

- **Move the Aquilex System, the other device, or both devices to a different location**
- **Increase distance between devices used**
- **Consult an electro-medical expert**

WARNING!
The Aquilex Fluid Control System should not be used directly next to other devices as this could result in malfunctions. The Aquilex Fluid Control System was tested for compliance with IEC 60601-1-2 as a stand alone system. Therefore, do not stack other devices (e.g. MyoSure® Control Unit) on the system or the Irrigation Pump Unit. In particular, do not place any other device than the AQL-100P on the trays of the AQL-100CBS. If usage in the manner described above is nevertheless required, this system and the other devices should be monitored to make sure they function properly.

WARNING!
If the Aquilex Fluid Control System is configured as part of a ME SYSTEM, the entire ME SYSTEM should be tested for compliance with IEC 60601-1-1, and any equipment used with the Aquilex Fluid Control System should be Type BF.

WARNING!
If the leakage current of the configured ME SYSTEM exceeds the limits of IEC 60601-1-1, install an appropriately rated UL 2601-1/IEC 60601-1 approved isolation transformer and retest the system.

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WARNING!
Always use the hooks of the bag scale to hang the fluid bags to ensure an accurate determination of the fluid deficit. In addition, leave the empty fluid bags hanging on the bag scale until the end of surgery.

3.1.2 Precautions

CAUTION!
Federal Law (only for U.S. market)
Federal law restricts this device to sale by or on the order of a physician.

CAUTION!
Indoor climate
Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

CAUTION!
When using the Aquilex System with tissue removal systems, e.g. MyoSure®, the combination of low set pressures and excessive vacuum pressures may result in a significant loss of intrauterine distension pressure which has the potential to affect the visibility of the surgical field. Conversely, when employing a high distension pressure, the deactivation of the tissue removal system can lead to pressure spikes that can exceed 150 mmHg.

CAUTION!
Do not use the covered power output at the rear of the irrigation pump unit.

CAUTION!
The system may only be connected with hysteroscopes designed for and featuring the technical specification permitting such a combined use. Any utilized hysteroscopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600.

CAUTION!
Check to ensure the available wall outlet voltage matches the data listed on the label attached to the back of the pump. Incorrect voltage can cause errors and malfunctions and may destroy the system.

CAUTION!
The device is transportable. The roller wheels of the Fluid Monitoring Unit (cart/scale) are used for positioning at the place of use. To transport the device, remove all fluid bags from the hooks and make sure there are no containers or only completely emptied containers on the cart/scale. Inflow and outflow tubes must be completely removed. Make sure the power supply line does not touch the ground and there are no other objects located on the Aquilex fluid control system. Always use the handle to move the system safely.

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CAUTION!
To avoid affecting the accuracy of the deficit calculation ensure that the first step of the container change is to disconnect tubing from the full containers. Remove full containers from the scale immediately after that.

Only for U.S. operators

Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 and IEC 16320-C13. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.

3.2 Description of the Aquilex Fluid Control System

Technical application scope of the system

The intrauterine pressure can be adjusted on the front of the pump. It can be pre-set to a range between 40 and 150 mmHg. The maximun inflow rate is 800 ml/min and is reduced automatically by the pump once the pre-set intrauterine pressure setting has been reached.

The system has been designed to provide both fluid and vacuum systems that maximize the performance of tissue removal systems, e.g. MyoSure®.

Suggested distension media

The Aquilex Fluid Control System can be used with hypotonic, electrolyte-free media (e.g., glycine 1.5% and sorbitol 3.0%) and isotonic, electrolyte containing media (e.g., saline 0.9% and Lactated Ringer's).

Pressure measuring and regulating

The system operates with a completely non-contact pressure measurement of the irrigation medium. The contact-free pressure measurement is achieved by integrating the pressure chamber into the tubing system. The pressure chamber transmits the irrigation fluid pressure to the electronics of the device via a pressure sensor. The pressure control circuit continuously compares the desired pre-set intrauterine pressure with the actual intrauterine pressure. The function of this algorithm is to maintain the pre-set intrauterine pressure. Check for possible leaks if the pre-set intrauterine pressure cannot be achieved.

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4 Initial System Set-up

Always check all parts and accessories of the system when performing initial set-up. If the system has obvious defects, contact Hologic Technical Support (Chapter 14, Warranty Information).

Initial system set-up

Place the system on a level surface and install in a dry environment. The ambient temperature and humidity must meet the requirements mentioned in Chapter 11, Technical Data.

WARNING!
Equipment should be positioned such that power cord can be easily disconnected.

CAUTION!
Indoor climate
Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

4.1 Preparing the System for Use

WARNING!
The Aquilex Fluid Control System should not be used directly next to other devices as this could result in malfunctions. The Aquilex Fluid Control System was tested for compliance with IEC 60601-1-2 as a stand alone system. Therefore, do not stack other devices (e.g. MyoSure® Control Unit) on the system or the Irrigation Pump Unit. In particular, do not place any other device than the AQL-100P on the trays of the AQL-100CBS. If usage in the manner described above is nevertheless required, this system and the other devices should be monitored to make sure they function properly.

Connection to the wall outlet

CAUTION!
Check to ensure the available wall outlet voltage matches the data listed on the label attached to the back of the system. Incorrect voltage can cause errors and malfunctions and may destroy the system.

Ensure the connection data and technical specifications of the power supply comply with DIN VDE or national requirements. The wall outlet power supply cord must be plugged into a properly installed safety wall plug (see DIN VDE 0107). Read the device label located in rear of pump to determine the operating voltage of the system.

Grounding contact

The power connection must be equipped with a grounding contact. Use the Aquilex power cord to establish a connection between the wall outlet and the power cord connection located in the rear of the system.

Only for U.S. operators

Use only a certified (UL-listed), removable power cord, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 or IEC 320/CEE22. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade outlet.

Potential equalization

Integrate the system into the potential equalization system as specified by local safety rules and regulations.

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en

Precautionary measures

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

This system is to be used only for the purposes described in the manual and has to be installed, set up, and operated in compliance with the EMC notes and instructions.

4.2 System Components


Fig. 4-1 System Components

1. Bag scale
2. Fluid bag
3. Irrigation Pump Unit
4. Trays
5. Container scale cable/connector
6. Container scale
7. Container holders
8. Container
9. Roller wheel base
10. Locking foot brake

The Aquilex Fluid Control System is divided into two separate boxes for shipping:

Box 1 contains:

- Irrigation Pump Unit
- Instructions for Use
- Power cord
- Aquilex system vacuum tube set (low and high vacuum)

Box 2 contains:

- Fluid Monitoring Unit (cart with scale)
- Container rings

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5 System Operation

Please make sure that the functional check according to chapter 6 has been performed prior to each device use.

5.1 Front of Irrigation Pump Unit


Fig. 5-1 Front of Irrigation Pump Unit

1. Pump display
2. Intrauterine pressure display
3. Fluid deficit limit display
4. Deficit meter
5. Deficit display
6. Inflow tube holder
7. Roller wheel
8. Pressure sensor
9. Reset deficit button (Zero)
10. Decrease deficit limit
11. Increase deficit limit
12. Decrease intrauterine pressure setting
13. Increase intrauterine pressure setting
14. Intrauterine pressure setting display
15. ON/OFF switch
16. Pause/Resume button
17. Prime button

Please familiarize yourself with the layout of the individual elements on the front of the pump.

5.2 Rear of Irrigation Pump Unit


Fig. 5-2 Rear of Irrigation Pump Unit

1. Connection for low vacuum (white)
2. Product label
3. Performance data of the device
4. Adjustment controller for high vacuum
5. Connection for high vacuum (green)
6. Fuse holder
7. Covered power output
8. Power connection for pump
9. Potential equalization connection
10. Service interface
11. Connection for scale
12. Suction openings

Please familiarize yourself with the layout of the individual elements at the rear of the pump.

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CAUTION!
Do not use the covered power output at the rear of the irrigation pump unit.

WARNING!
Additional equipment
Additional equipment connected to medical electrical devices must be demonstrated to be compliant with their respective IEC or ISO standards (IEC 60601-1, IEC 60950 or IEC 62368 for data processing equipment).

5.3 Fluid Monitoring Unit Set-up


Fig. 5-3 Fluid Monitoring Unit (cart with scale)

1. Bag scale
2. Fluid bags
3. Bag deflector
4. Container scale cable/connector
5. Container scale
6. Holder for upper container (Serres, Medela)
7. Holder for lower container (Abbott, Bemis®, Medi-Vac®, DeRoyal®)
8. Roller wheel base
9. Container

The fluid monitoring unit (cart/scale) consists of a weighing unit for fluid bags (1), a weighing unit for fluid containers (5), and a roller wheel base (8).

1. Remove the cart/scale from the cardboard shipping box.
2. Remove the pump and the power cords from the first cardboard box.
3. Loosen the handwheel (6) (Fig. 5-4, page 17) and pull bag scale upwards to stop. The screw (7) (Fig. 5-4, page 17) must be inserted into the provided opening. Secure the bag scale with the handwheel.
4. Depending on the type of container used, attach the container rings (included in the second box) to the upper (6) or lower (7) container holders (Fig. 5-3, page 16).
5. Guide the power cord through the holes provided for this purpose and connect to pump (2) (Fig. 5-4, page 17) and plug into grounded, shockproof safety wall socket.
6. Attach the scale to the pump by connecting the connector of the bag scale with the connector of the container scale (4) (Fig. 5-4, page 17) and fix the connected cables below the lower pump tray by means of the provided wire clips. Connect the second connector of the bag scale to the port on the back of the pump (5) (Fig. 5-4, page 17).

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
Fig. 5-4 Scale and pump connection

1. Bag deflector
2. Power cord/pump connection
3. Covered power output
4. Plug connection, bag scale with container scale
5. Bag scale connector, connection to pump
6. Handwheel
7. Screw

NOTE!
The pump's software automatically detects whether the container scale and bag scale are connected or only the container scale.

WARNING!
Scale error
Ensure that nothing weighs down the scale during system start-up. Doing so may result in an inaccurate deficit value.

WARNING!
Fluid deficit
The fluid left in the patient must be monitored. The deficit is the total amount of fluid left in the patient or unaccounted for otherwise. Take notice of the measurement tolerance of the system (see Chapter 11, Technical Data). Estimating the fluid volume remaining in the patient is the physician's responsibility.

NOTE!
The greater the consumption of irrigation fluid, the greater the deviation between the actual and the displayed deficit (see "Technical Data," deficit accuracy: ± 6%).

WARNING!
Serum sodium concentration
It is also necessary to monitor the concentration of sodium in the blood of the patient to prevent electrolyte disturbances. Monitoring of the concentration of sodium in the blood must be performed by the physician and is not performed or supported by the system.

Precise balancing

Try to collect all the fluid running out of the uterine cavity during the procedure to achieve the most exact deficit value possible.

Scale capacity

The container scale can be loaded with a weight of up to 25 kg (55 lbs). The maximum load of the bag scale is 12 kg (27 lbs). If the scales are loaded beyond these

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






limits, the message **Scale Overload/Check Scale** is triggered and displayed. Three audible warning beeps are emitted as well (See Chapter 10, Error and Warning Messages.)

CAUTION!
Make sure the containers and fluid bags hang freely, are not resting on something, and do not touch other objects except the bag deflectors. Failure to follow these instructions means the deficit cannot be calculated correctly.

NOTE!
Connect the scale to the pump before turning the system on to ensure the system recognizes the scale.

5.3.1 Setting of the Container Scale

The container scale can be used with containers from different manufacturers.

Bemis®3 liters		DeRoyal® Crystaline™ 2.1 l	
Abbott 2 liters		Serres 2 & 3 liters	
Medi-Vac® 3 liters		Medela 3 liters	
Medi-Vac Flex Advantage 3000 cc			

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NOTE!
Ensure containers are positioned properly in the respective holders.

NOTE!
Only use containers with overflow protection.

5.3.2 Connecting the Vacuum Tube

CAUTION!
When using the Aquilex System with tissue removal systems, e.g. MyoSure®, the combination of low set pressures and excessive vacuum pressures may result in a significant loss of intrauterine distension pressure which has the potential to affect the visibility of the surgical field. Conversely, when employing a high distension pressure, the deactivation of the tissue removal system can lead to pressure spikes that can exceed 150 mmHg.

Connect vacuum tube with hygiene filter to suction containers. The vacuum tube with hygiene filter must be replaced when dirty and after 30 days at the latest. The vacuum tube with hygiene filter should not be cleaned.

- Connection for low vacuum (white)

▶ Connect vacuum tube with white connector to low vacuum port (white) (1) Fig. 5-5. This vacuum pump has a fixed vacuum pressure (~ 225 mmHg).

▶ Use the connecting tube ((5) Fig. 5-5) when two containers are serially connected to the same vacuum port.


Fig. 5-5 Low vacuum tube

(1) Connection for low vacuum (white)
(2) Hygiene filter
(3) Vacuum tube
(4) container
(5) Connecting tube

- Connection for high vacuum (green)


▶ Connect vacuum tube set with the green connectors to the high vacuum port (green) (8) in Fig. 5-6. This vacuum can be adjusted to a maximum 500 mmHg using adjustment dial.

▶ Use the connecting tube ((12) Fig. 5-6) when two containers are serially connected to the same vacuum port.

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
Fig. 5-6 High vacuum tube

- (6) Hygiene filter
- (7) Vacuum tube (green connectors)
- (8) Connection for high vacuum (green)
- (9) Container
- (10) Tissue trap
- (11) Adjustment dial
- (12) Connecting tube

5.4 Turning On the Aquilex System

1. Press the ON/OFF switch. The displays and indicators light up and system turns on.
2. The system now performs a device self-test.

 The device self-test is used by the system to check whether the container and the bag scale are connected, among others. If the bag scale is not connected, the message "System OK. Bag Scale Not Connected" is displayed. In this case, please check the appropriate data connection of the scales with the pump according to chapter 5.3.
3. If a tube set is in the inflow tube holder when the pump is switched on, the **pump display** (Fig. 5-1, Front of Irrigation Pump Unit (1)) shows the message **Remove Tube Set**. The device self-test resumes once the tube set is removed from the roller wheel.

 If the device self-test is unsuccessful, the corresponding error messages are displayed (see Chapter 10, Error and Warning Messages).

The system has successfully completed the device self-test when a single audible beep is heard. The message **System OK** is displayed for 5 secs followed by the message **Insert Tube Set**.

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WARNING!
Do not use this system if a defect is suspected or detected during the function check. This also applies to obvious defects, especially defects and damage to the power plug and power cord.

5.5 Hanging the Fluid Bags


Fig. 5-7 Fluid bag suspension

① Bag scale
② Fluid bag on bag hook
③ Bag deflector

WARNING!
When performing monopolar hysteroscopic electrosurgery, the distension medium must be electrically non-conductive. Examples include glycine, sorbitol and mannitol. Isotonic saline irrigation fluids may only be used when performing bipolar electrosurgical resective procedures.

Hang one or two fluid bags filled with distension fluid appropriate for procedure.

WARNING!
The system is only intended for use with flexible fluid containers. Do not use glass containers as they might break. With rigid containers, fluid cannot flow quickly enough due to the vacuum being generated inside of the containers. Risk of implosion with rigid containers.

5.6 Using Tube Sets

The Aquilex Fluid Control System is designed for use with sterile disposable inflow and outflow tube sets.

Tube set recognition

Each inflow tube set is equipped with tube set recognition technology. An RFID transponder detects the type of tube, whether it has been used, and its reliability automatically. The pump display indicates this information. This eliminates accidental reuse of tube sets on more than one patient (see Chapter 5.7, Tube Overview).

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WARNING!
Visual inspection of the tube set
Before the operation, perform a visual inspection of the tube set and its packaging. Damaged tube sets or tube sets from damaged packagings may not be used.

WARNING!
Reprocessing of sterile disposable products
Reuse of inflow or outflow tube can cause an infection hazard for patients and/or users as well as impair of product functionality. Contamination and/or impaired functionality of the system can cause risk of injury, illness, or death. Do not re-process or reuse single-use inflow or outflow tube sets.

NOTE!
Comply with national disposal and hygiene rules when disposing of tubes, collected fluid, and the containers or containers.

5.7 Tube Overview

Three different tube sets are necessary to operate the system. The following table lists each type of tube set and its application.

Article number	Description
AQL-110	Aquilex Fluid Control System inflow tube set, disposable, sterilized using ethylene oxide
AQL-111	Aquilex Fluid Control System outflow tube set, disposable, sterilized using ethylene oxide
AQL-112	Complete tube set (inflow and outflow) disposable, sterilized using ethylene oxide
AQL-114	Aquilex Fluid Control System high and low vacuum tube set: re-usable, non-sterile

Table 5-1

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5.8 Connecting the Outflow Tube Set

CAUTION!
When using the Aquilex System with tissue removal systems, e.g. MyoSure®, the combination of low set pressures and excessive vacuum pressures may result in a significant loss of intrauterine distension pressure which has the potential to affect the visibility of the surgical field. Conversely, when employing a high distension pressure, the deactivation of the tissue removal system can lead to pressure spikes that can exceed 150 mmHg.



Fig. 5-8 Outflow tube set

1. To low vacuum port (white)
2. Container
3. Connecting tube
4. Patient port
5. Outflow tube set
6. Drape
7. Removable outflow channel or hysteroscope outflow sheath stopcock


Using low vacuum configuration of Fig. 5-8, connect outflow tube set (Y-tube) to patient port (4) of second container. Yellow flexible connector attaches to drape (6). Yellow luer fitting connects to stopcock (7) of removable outflow channel or hysteroscope outflow valve.

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
Fig. 5-9 Port for tissue removal systems

① To high vacuum port (green)
② Container
③ Connecting tube
④ Specimen tissue port
⑤ Tissue trap
⑥ Vacuum tube of tissue removal handpiece (yellow)
⑦ Tissue removal handpiece

5.8.1 Connecting Outflow Tube of Tissue Removal Handpiece (e.g. MyoSure®)

If intrauterine pathology is identified, the outflow tube of a tissue removal handpiece ⑥ is connected to the tissue trap ⑤ located in the second container.

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5.9 Inserting the Inflow Tube Set


Fig. 5-10 Tube set elements

1. Protective caps
2. Bag spikes
3. Tube clamps
4. Y-connector
5. Inflow section
6. Roller wheel section
7. Pressure chamber with membrane and RFID transponder
8. Hysteroscope section
9. Luer lock connector (blue)
10. Roller wheel connector

(See Fig. 5-10, Tube set elements) The inflow tube set consists of three tube sections, a Y-connector (4) and two bag spikes (2). The three tube sections are: Roller wheel section (6), inflow section (5), and hysteroscope section (8). The bag spikes (2) are used to connect the tube sections to the fluid bags.

The luer lock connector (9) connects the hysteroscope tube with the hysteroscope.


Fig. 5-11 Inserting the tube set

1. Bag spikes
2. Fluid bags
3. Bag clamps
5. Inflow tube
6. Roller wheel tube
7. Pressure chamber with membrane and RFID transponder
8. Hysteroscope tube
9. Luer lock connector (blue)

1. Inflow tube set - To be carried out by non-sterile technician:

Open outer packaging

- ▶ Open outer packaging of the inflow tube set.
- ▶ A sterile nurse then removes the inner tube set package and opens it.

2. To be carried out by sterile technician:

Connect to hysteroscope

- ▶ Keep the blue luer lock connector (9) in the sterile area and hand the tube end with the bag spikes (1) to the non-sterile technician.

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Inserting the tube set

- ▶ Connect the blue luer lock connector (9) with the hysteroscope inflow stopcock. Open stopcock.

3. To be carried out by non-sterile technician:


- ▶ Ensure system is turned on.
- ▶ Close the clamps (3) on the inflow tube below the bag spikes (1).
- ▶ Insert the inflow tube set into the inflow tube holder.
 Insertion of the roller wheel tube is depicted in Fig. 5-12.
- ▶ Carefully insert the pressure chamber (7) into the lower notch of the inflow tube holder (12) until you feel resistance. Align pressure chamber and inflow tube holder using arrows (see Fig. 5-12).
- ▶ When inserting the roller wheel tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber (7) only if chamber is not pressurized.
- ▶ Place the roller wheel tube (6) around the roller wheel (11).

Connect the fluid bags

- ▶ When connecting or removing the tube to or from the irrigation fluid bags, always grasp the bag spike at the provided handle. Observe aseptic technique when inserting the spike(s) into the fluid bag(s). The surgeon must select a distension fluid suitable for the type of procedure.

Fig. 5-12 Attach roller wheel tube

(5) Inflow tube
(6) Roller wheel tube
(7) Pressure chamber
(8) Hysteroscope tube
(11) Roller wheel
(12) Holder for inflow tube
(13) Alignment arrows

5.10 Presetting the Intrauterine Pressure

Intrauterine pressure setting

The intrauterine pressure setting can be adjusted while the system is in operation. Use the ▲ and ▼ buttons (Fig. 5-1, Front of Irrigation Pump Unit). The pressure setting can be adjusted to between 40 to 150 mmHg in steps of 5 mmHg. The intrauterine pressure is shown on the intrauterine pressure display (2).

Safety threshold

If when scrolling with the ▲ button (Fig. 5-1, Front of Irrigation Pump Unit) the safety threshold of 100 mmHg is reached, an audible warning beep is emitted. Release the ▲ button for one second and scroll again to set higher values up to 150 mmHg.

CAUTION!
If the intrauterine pressure does not react to an increase in the pressure setting during the procedure, a perforation of the uterine cavity might be the cause. This results in an increased risk of intravasation. Examine the uterine cavity for injuries.

5.11 Deficit Limit Setting

Setting the Deficit Limit

The deficit limit can be adjusted while the system is in operation. Use the ▲ and ▼ buttons (see Fig. 5-1). The deficit limit can be adjusted to between 600 and 2500 ml in increments of 100 ml. The deficit limit is shown on the deficit limit display (3). The deficit meter is designed to help the operator track the deficit volume. The color of the deficit meter changes as the deficit limit is approached. The deficit limit set by the operator is marked with a red LED on the top of the deficit meter. If during surgery the actual deficit rises, the LEDs will light up sequentially,

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representing the actual deficit volume until the deficit limit is reached (see section **Deficit Limit** in Chapter 7, Safety Functions).

5.12 Using the Pump during Surgery

- ▶ Open bag clamps on the fluid bags (③ Fig. 5-11).
- ▶ Fully open hysteroscope inflow stopcock.
- ▶ If drainage stopcock is available: Fully close drainage stopcock.
- ▶ Keep the hysteroscope at the height of the patient and above the drape to ensure the fluid can be collected. Do not insert the hysteroscope into the uterus at this time.
- ▶ Press the button **Prime** (⑰ Fig. 5-1).
- ▶ Pump will run for approximately 20 seconds to purge air from tubes and run the automatic lumen calibration function.
- ▶ Pump will display **Calibration Running.**

Automatic lumen calibration

The pump is equipped with an automatic lumen calibration function. The system determines the flow resistance of the hysteroscope. This resistance is used to calculate the pump pressure necessary to maintain the pre-set intrauterine pressure. To overcome this resistance, the pump allows a pressure of up to 80 mmHg during calibration. This will be shown in the display of the actual intrauterine pressure. If calibration fails due to high resistance, calibration is repeated with a permissible pressure of up to 150 mmHg. If calibration is then still not completed, the pump displays **Prime Fail - Open Stopcock Clamps**.

The automatic lumen calibration starts once the **Prime** button is pressed.

- ▶ Three beeps are heard once the automatic lumen calibration is finished. The pump display will show **Prime Successful Close Stopcock** for 5 seconds followed by **System Operating.**
- ▶ Close hysteroscope inflow stopcock to stop fluid flow. Once all fluid has been removed from the drape, zero the deficit display.
- ▶ Check to see if fluid has leaked in the area of the pressure chamber. If you find leaked irrigation fluid in the area of the pump, change the tube set and retry the automatic lumen calibration.

NOTE!
The calibration must be performed outside of the patient to ensure a correct lumen calibration and deficit calculation.

NOTE!
The pump continues to operate after automatic lumen calibration is complete. The pump should be stopped by closing hysteroscope inflow stopcock.

CAUTION!
A new lumen calibration must be carried out with each change of the hysteroscope during surgery (5.15, Changing Instrument during Surgery).

System Operation

- ▶ Open stopcock and guide the hysteroscope with fluid flowing into the uterus.
- ▶ Adjust intrauterine pressure setting as necessary to obtain adequate distension and visualization.

Completing System Operation

- ▶ After the medical procedure, if system operation is complete close the hysteroscope inflow stopcock.
- ▶ Wait until the entire fluid volume from the drape and the tube set has been collected in both containers.
- ▶ Press the **Pause/Resume** button.
- ▶ Note the deficit volume indicated on the deficit display. This is the total fluid volume that was absorbed by the patient.

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WARNING!
System error: Do not use the Aquilex System if a defect is suspected or detected during the function check. This also applies to any obvious defects, especially defects on the power connector or plug and power cord.

WARNING!
If the message "Check Scale Connection" appears, the deficit must be calculated manually. The pump keeps displaying the last known deficit value determined prior to the failure of the scale connection.

NOTE!
It is possible to change the containers and bags during surgery without losing the previously measured deficit.

5.13 Changing Bags during Surgery

Bag change during surgery

The device system detects automatically the replacement of a bag. Brief fluctuations in the deficit calculation (< 10 s) may occur when a bag is replaced. The replacement of a bag is indicated with the message **Bag Change, Please Proceed**.

- ▶ Close the clamp of the empty bag.
- ▶ Hang a new fluid bag to a hook of the bag scale.
- ▶ Reconnect the new fluid bag with the inflow tube set.

WARNING!
Touching the bags and the hooks as well as vibrations of the balancing system should be avoided during surgery to prevent false detection of the bag change and not negatively affect the accuracy of the deficit calculation.

WARNING!
The empty bags should remain on the hooks so that the measurement accuracy is not diminished.

WARNING!
Bags should be changed quickly to avoid affecting the accuracy of the deficit calculation.

5.14 Changing Container during Surgery

Container change during surgery

The device system detects automatically the replacement of a container. The pump will stop immediately and the deficit display is locked to insure an accurate deficit count is maintained. Brief fluctuations in the deficit calculation (< 10 s) may occur when a container is replaced. The replacement of a container is indicated with the message **Container Change, Press Resume**.

- ▶ Disconnect tubing from full containers.
- ▶ Immediately remove full containers from scale.
- ▶ Install new containers.
- ▶ Reconnect tubing to new containers.
- ▶ Press **Pause/Resume** button to resume procedure.

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WARNING!
A container change during surgery is only allowed, if the container holds at least 0,5 liters of fluid. Otherwise, the deficit value may be falsified. In this case, the manufacturer recommends manual deficit calculation.

WARNING!
Touching the containers and their holders as well as vibrations of the balancing system should be avoided during surgery to prevent false detection of the container change and not negatively affect the accuracy of the deficit calculation.

WARNING!
Containers should be changed quickly to avoid affecting the accuracy of the deficit calculation.

CAUTION!
To avoid affecting the accuracy of the deficit calculation ensure that the first step of the container change is to disconnect tubing from the full containers. Remove full containers from the scale immediately after that.

5.15 Changing Instrument during Surgery

Instrument change during surgery

- ▶ Pause the pump by pressing **Pause/Resume** button.
- ▶ Press the **Prime** button for 2 seconds.
- ▶ Change the instrument.
- ▶ Fully open hysteroscope inflow stopcock.
- ▶ Keep the hysteroscope at the height of the patient and above the drape to ensure the fluid can be collected. Do not insert the hysteroscope into the uterus at this time.
- ▶ Press the **Prime** (17) button (Fig. 5-1).
- ▶ Pump is running to carry out the automatic lumen calibration. The pump display depicts **Calibration Running**.
- ▶ Three beeps are heard once the automatic lumen calibration is finished.
- ▶ The pump display depicts **Prime Successful Close Stopcock** for 5 seconds, followed by **System Operating**.
- ▶ Close the stopcock for the hysteroscope inflow to stop the inflow.

5.16 Total Inflow Volume Displayed

Total inflow volume displayed

If a manual check of the fluid deficit is desired, the total fluid volume supplied from the fluid bags can be obtained by simultaneously pressing and holding both the up and down arrows ((10) and (11) Fig. 5-1). The number in the Deficit Display is the total inflow fluid volume in ml. Once one or both of these buttons is released, the Deficit Display will return to the Fluid Deficit value.

5.17 Turning System Off

Turning Off

Press the **ON/OFF** switch to turn pump off. The displays and indicators are no longer illuminated.

WARNING!
Pressing the ON/OFF switch does not disconnect the system from the wall power outlet. This requires pulling the power cord located in the rear of the system.

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6 Functional Check

WARNING!
Functional test
The functional test must be performed prior to each device use.

1. Perform a visual check of the devices. Do not use the system in case of obvious damage.
2. Check the rollers of the roller wheel to make sure they move easily and smoothly.
3. Turn on the device, check whether power switch and indicators and displays light up accordingly.
4. The device self-test must run successfully; there are no displayed error messages (5.4, Turning On the Aquilex System).
5. The bags with the irrigation fluid must hang freely and may not touch the scale.
6. Check that all tube connections (vacuum/inflow/outflow) were made correctly and are intact.
7. Check to make sure all tube connections are free of mechanical stresses and are routed without snagging. The tube connections may not touch the scale. Non-observance can lead to a distortion of the deficit calculation.
8. The automatic lumen calibration has been carried out successfully; there are no displayed error messages (5.12, Using the Pump during Surgery).
9. Check for leaking irrigation fluid in the area of the pressure chamber.

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7 Safety Functions

The electronic components continuously monitor the proper function of the system. System malfunctions are indicated with audible warning beeps, error messages, and/or the blocking of system functions. A table listing a summary of possible error and warning messages is provided in Chapter 10, Error and Warning Messages.

Intrauterine pressure 10 mmHg above preset intrauterine pressure setting

If the intrauterine pressure exceeds the intrauterine pressure setting by 10 mmHg for longer than 5 seconds, the pressure reduction function is activated. The roller wheel will move forward or backward a few times during the pressure reduction process. If the pressure cannot be reduced, the message **Overpressure/ Open Stopcock** is displayed and three warning beeps are emitted.

Intrauterine pressure > 150 mmHg

The message **Maximum Pressure** is displayed and the pump will emit 3 warning beeps once the intrauterine pressure exceeds 150 mmHg. The maximum permissible pressure has now been reached.

Intrauterine pressure > 200 mmHg

If the intrauterine pressure exceeds 200 mmHg for longer than 5 seconds, the roller wheel stops and the message **Overpressure/Check Stopcock** is displayed. Three short continuous audible warning beeps are emitted until the pressure is reduced. Once the intrauterine pressure falls below 200 mmHg, the audible warning beeps stop and the pump wheel resumes turning automatically.

Check Tube Set Installation

If the inflow tube set is not inserted properly into the roller wheel, pressing the **Prime** button results in a short audible beep and **Check Tube Set Installation** is displayed. The roller wheel does not start to turn.

Pressure measuring system malfunctions

If a malfunction is detected in the pressure measurement electronics, **Sensor Error** is displayed and five short warning beeps are emitted. The roller wheel stops turning.

Scale overload

If the maximum permissible weight of the scale is exceeded (bag scale or container scale, a continuous audible warning beep is emitted and **Scale Overloaded Check Scale** is displayed. The warning stops once the excess weight is removed from the scale.

Loading/unloading scale while in operation

The following messages are displayed when a container is removed from the scale or a new fluid bag is added during operation of the system: **Container Change, Press Resume or Bag Change, Please Proceed**.

Pressure setting at restart

If the last pre-set intrauterine pressure setting is greater than 80 mmHg, this value is reset to the default value of 80 mmHg.

Deficit limit

Each additional deficit increase by 100 ml above the selected deficit threshold triggers 3 warning beeps that are repeated while the pump continues to operate. The message **Deficit Limit Exceeded** is displayed.

Deficit rate >300 ml/min

If the deficit rate exceeds 300 ml/min, three audible warning beeps are emitted and the message **High Fluid Loss Check Leakage** is displayed. If no obvious source of high fluid loss can be identified, an assessment of potential cervical or uterine perforation should be made.

Serious system defect

Five short audible warning beeps are emitted and the message **Motor Error** is displayed. It is possible that errors can occur in the start up sequence prior to the enabling of the Pump Display. In this situation, the Pump Display will remain blank.

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Care and maintenance

8 Care and Maintenance

The service and maintenance of the system and its accessories has to be carried out as per instructions to ensure the safe operation of the system. For the protection of the patient and the operating team, check that the system is complete and functional before each use.

Special care is necessary when servicing, maintaining, and storing the system and its accessories to maintain the functionality of the equipment and any attached devices.

8.1 Cleaning the System

1. Use the **ON/OFF** switch to turn off the system.
2. Remove the power cord.
3. Wipe the surface of the system with a soft cloth moistened with a disinfectant (for example, Meliseptol® rapid). The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the system.

WARNING!
Cleaning the system / Sterilization not allowed.
The pump and the cart/scale can be disinfected by wiping off the outer surfaces.
Do not sterilize the pump and the cart/scale.

8.2 Maintenance Carried out by Authorized Service Technician

Two-year maintenance interval

It is recommended that an authorized service technician inspects and services the system at appropriate intervals to ensure safety and functionality. The **minimum** service interval is two years, depending on frequency and duration of use. If this interval is not maintained, the manufacturer does not assume any liability for the functional safety of the system. A sticker located on the rear panel of the system contains the latest date for the next service or maintenance check.

NOTE!
Service and maintenance work may not be carried out during surgery.

Certification

Ask the service technician for a certificate after he or she has inspected the system or performed any service tasks. This certificate must list:

- type and scope of service,
- date of service,
- name of company performing service
- as well as signature.

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8.3 Replacing of the Fuse

CAUTION!
Before replacing the fuse, check the values of the fuse to be inserted according to Chapter 11, Technical Data, page 42.

The fuse may be defective and is in need of replacement if:

- one or more of the pump displays does not light up,
- the system does not function.

Check whether

- the power cord is properly connected to the power cord connection (Figure 5.2) and to a grounded safety wall outlet,
- the wall outlet has power.

WARNING!
Unplug the power cord from the system before checking the fuse.

The system does **not** have to be opened to replace the fuse.

1. Turn system off.
2. Disconnect system from wall power outlet.
3. The fuse holder is located on the back of the pump, next to the male connection.
4. Remove both fuse holders as depicted in Fig. 8-1, using small flathead screwdriver.
5. Pull out the fuse holders.
6. Check the fuses.
7. Insert new fuses. Use only the specified type of fuse (see Chapter 11, Technical Data).
8. Insert the fuse holders.
9. Reconnect the power cord and connect the pump to the wall outlet.


Fig. 8-1 Opening the fuse holder

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Manufacturer's specification

9 Annual Inspection

The manufacturer stipulates that qualified personnel or biomedical technicians must regularly test the system to assess its functionality and technical safety. These inspections must be carried out annually. Regular inspections will assist in early detection of possible malfunctions. This helps maintain the system and increases its safety and service life.

Inspection tests

The following tests are designed specifically for trained personnel or a biomedical technicians. System operation, functionality, and serviceability are easily checked. Each test conducted must be documented by signing and dating the test log.

WARNING!
If the specified parameters and tolerances are exceeded, the system must be returned to Hologic for evaluation.

9.1 Safety Test

1. Perform a visual inspection. Ensure:
 - the fuse corresponds with the specifications indicated by the manufacturer,
 - labels and stickers on system are legible,
 - the mechanical condition of the device allows for its safe use,
 - the device is clean to ensure proper and safe functionality.
2. Perform the measurement of the ground leakage current (max. 500 µA) and contact current (max. 100 µA in normal state and max. 500 µA on single fault condition) according to IEC 60601-1 / EN 60601-1.
3. Measure protective conductor resistance according to IEC 60601-1 / EN 60601-1. The protective conductor resistance is to be measured with a connected power cord; however, the cord may not be connected to the mains power supply. The max. value is 0.2 Ω.

As an alternative, perform safety test according to DIN EN 62353. Measuring the protective conductor resistance should be carried out only according to IEC 60601-1 / EN 60601-1 using a current of 25 A AC.

9.2 Basic Function Tests

The basic function tests check displays, buttons, and overall system performance. You need the following for this test:

- Aquilex inflow tube set
- Fluid bags
- Measuring cup with marked scale of 1 liter
- Stopwatch
- Precision weights (e.g., Ohaus 1 kg 49016-11 or 41000-00 or equivalent)

NOTE!
If the device does not work as described and the test fails the device must be sent to the service.

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9.2.1 Scale Test

1. Switch system on.
2. Once the display depicts **Insert Tube Set**, press **Pause/Resume** and **Zero** at the same time.
3. The pump display depicts the message **Scale Test**.
 The first option is the test of the container scale.
4. Place a precision weight on the container scale (500 g - 2000 g)
5. The display of the fluid deficit limit depicts the weight.
6. The tolerance range is ± 20 g.
7. If the displayed value is outside the tolerance range, the scale must be calibrated by a service technician.
8. Remove weight from container scale.
9. Press the Deficit Level UP button to start the test of the bag scale.
10. Place a precision weight on the bag scale (500 g - 2000 g).
11. The difference weight is depicted in the deficit display.
12. The tolerance range is ± 20 g.
13. If the displayed value is outside the tolerance range, the scale must be calibrated by a service technician.
14. Remove weight from bag scale.
 (Switch to the container scale test by pressing the Deficit Level Down button)
15. Press the **Pause/Resume** button to end the test.

Enter the results into the test log in section 16.1. The test is successfully completed when the results are within the permissible tolerance limit.

9.2.2 Flow Rate Test


Fig. 9-1 Flow Rate Test

The test setup is depicted in Fig. 9-1, Flow Rate Test.

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1. Switch system on. (See Chapter 5.4, Turning On the Aquilex System)
2. Insert tube set into pump and close clamps on the bags.
3. Hang the fluid bags onto the fluid bag hooks.
4. Insert tap spikes into the fluid bags and open clamps on the fluid bags.
5. Insert hysteroscope tube into the measuring cup.
6. Set intrauterine nominal pressure to 150 mmHg.
7. Press the **Prime** button
8. The roller wheel starts to rotate in order to remove the air from the tubes and to carry out the automatic lumen calibration.
9. When the automatic lumen calibration is completed, press **Pause/Resume** button (about 20 seconds).
10. Empty the measuring cup.
11. Place hysteroscope tube back into the measuring cup.
12. Press the **Pause/Resume** button.
13. Press **Pause/Resume** button after one minute. The measuring cup should contain approx. 800 ml +/- 60 ml of fluid.

Enter the results into the test log in section 16.1. The test is successfully completed when the results are within the permissible tolerance limit.

9.2.3 Pressure Measuring Test

The test setup is depicted in Fig. 9-2.


Fig. 9-2 Set-up of pressure measuring test

h Height of the water line

The pressure test checks the pressure chamber, the pressure sensor, and the accuracy of the pressure measurement, to ensure that all elements function correctly. For this test, an inflow tube set and a **water-filled** container is required. The height of the water column (hydrostatic pressure) is used to test the pressure transducer.

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1. Hang the inflow end of the tube with the tap spikes for the fluid bags into a container filled with water.
2. Fill the tube end completely with water by running the pump with the **Prime** button. Let the pump run until the calibration is complete. Press the **Pause/ Resume** button to stop the roller wheel. The intrauterine actual pressure display depicts 0 mmHg.
3. Seal the end of the hysteroscope tube (with finger on the tip of the luer connector).
4. Hold the water level at the end of the hysteroscope tube (h) 30 cm above the pressure chamber. The water column exerts a hydrostatic pressure on the pressure transducer.
5. Remove your finger from the end of the hysteroscope tube.
6. The intrauterine actual pressure display should depict 20 mmHg (±5 mmHg).
7. Change the water column height by changing the height of the water filled end of the tube set. The value depicted on the intrauterine actual pressure display should change accordingly.

Enter the results into the test log in section 16.1. The test is successfully completed when the results are within the permissible tolerance limit.

9.2.4 Fluid Deficit Measurement Test


The test setup is depicted in Fig. 9-3. It is critical that the collection container is positioned **on** the scale as shown in Fig. 9-3.

1. Carry out the scale test, chapter 9.2.1
2. Set the fluid deficit display to zero by pressing the **Zero** button (see Fig. 5-1, ⑨).
3. Press the **Pause/Resume** button.
4. Allow system to run for 1 minute.
5. Press the **Pause/Resume** button.
6. The displayed fluid deficit should be 0 ml. The permissible tolerance is ±50 ml.

Enter the results into the test log in section 16.1. The test is successfully completed when the results are within the permissible tolerance limit.

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Fig. 9-3 Test Setup for Fluid Deficit Measurement

9.2.5 Testing the Vacuum Pump


This test is not designed as a performance test for the measurement of negative pressure. The test only indicates whether or not the vacuum pump is operational.

1. Check to be sure at least one port is opened into the container.
2. Press the **Pause/Resume** button.
3. Place a finger lightly over the suction openings (Fig. 9-4, ①) and check whether you can feel a sucking air flow.

Enter the results into the test log in section 16.1. The test is successfully completed when a stream of air can be demonstrated.

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Fig. 9-4 Vacuum pump exhaust ports

9.3 Determine the Software Version

The software version of the pump may be required for further analysis.

Follow these steps to determine the software version:

1. Turn device on and wait for the device self-test.
2. Press the Pause/Resume button for at least 2 seconds.
3. The software version of the pump is depicted in the "Deficit Limit" display and consists of a 5-digit numerical sequence.

Enter the determined software version into the test log in Section 16.1.

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10 Error and Warning Messages

The messages appear in the pump display and warning beeps are heard. Beeps indicating a warning or operating message are sounded several times.

Message in the pump display	Warning beeps	Method
Check Tube Set Installation	1 warning beep	Remove and re-insert tube set. If the message reappears, insert a new tube set.
Tube Set Over Usage Limits	1 warning beep	Tube set detection indicates that the tube set has been used already. Insert new tube set.
Check Flow Path, Stopcock, Clamps	3 warning beeps	Flow path is blocked. Check that bag clamps and hysteroscope stopcock are open. Make sure the tube set is not blocked.
Incorrect Tube Set	1 warning beep	Replace tube set. The tube set does not match the type approved for the Aquilex System.
Pump Paused, Press Resume	1 warning beep	Pause/Resume button has been activated. Press Pause/Resume button again to resume surgery.
Overpressure Open Stopcock	3 warning beeps	Most commonly triggered when hysteroscope stopcock is closed while pump is operating at peak flow rate. Open hysteroscope stopcock or remove other closure to relieve pressure.
Overpressure Check Stopcock	3 beeps, constantly repeating until pressure is reduced	Pressure has exceeded the 200 mmHg safety limit and must be reduced. The most common cause is a closed hysteroscope stopcock while the pump operates at the highest flow rate. Open inflow stopcock on hysteroscope or remove blockage clamping off the inflow tube set.
Maximum Pressure	No beep	The message "Maximum Pressure" is depicted when the intrauterine pressure is greater than 150 mmHg.
Deficit Limit Reached	3 warning beeps	Actions are at the discretion of the physician.
Deficit Limit Exceeded	3 warning beeps	Actions are at the discretion of the physician. If necessary, perform manual deficit control.
Pressure Threshold	1 warning beep	For a hysteroscopy, pressures above 100 mmHg are usually not necessary. Careful monitoring of the fluid deficit is recommended.
Check Scale Connection	3 warning beeps	Check scale connection. Reconnect scale, restart device. If the message reappears, contact Hologic.
Remove Tube Set for System Check	1 warning beep	Make sure all tube sets are removed from the roller wheel during the system test. Remove tube set and wait until you hear a warning beep and the message "Insert Tube Set" appears.
Prime Fail-Open Stopcock, Clamps	3 warning beeps	Check bag clamp(s) and hysteroscope inflow stopcock are open. Press the "Prime" button to restart.
High Fluid Loss Check Leakage	3 warning beeps	Actions are at the discretion of the physician. If necessary, perform manual deficit control.

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Message in the pump display	Warning beeps	Method
Scale Overload/Check Scale	3 warning beeps	Weight on the **container scale** exceeds 25 kg (55 lbs) or weight on the **bag scale** exceeds 12 kg (26.5 lbs). The weight must be reduced. System function will resume when the excess weight has been removed.
Communication Error	5 warning beeps	Contact Hologic Technical Support.
Calibration Error	5 warning beeps	Contact Hologic Technical Support.
Sensor Error	5 warning beeps	Contact Hologic Technical Support.
Motor Error	5 warning beeps	Contact Hologic Technical Support.
Low Vac Failed Use Alternative	3 warning beeps	An alternative low vacuum source must be used in order to continue the procedure. Contact Hologic Technical Support.
High Vac Failed Use Alternative	3 warning beeps	An alternative high vacuum source must be used in order to continue the procedure. Contact Hologic Technical Support.
Vac Systems Out Use Alternative	3 warning beeps	An alternative vacuum source must be used in order to continue the procedure. Contact Hologic Technical Support.

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en

11 Technical Data

Type designation	Aquilex Fluid Control System (REF: AQL-100S) consisting of: Irrigation Pump Unit (REF: AQL-100P) Fluid Monitoring Unit (REF: AQL-100CBS)	
Manufacturer Information	W.O.M. WORLD OF MEDICINE GmbH Salzufer 8, 10587 Berlin	
Software version	See user manual/operator instructions to determine software version (9.3, Determine the Software Version)	
Mains voltage range [V]	100-240 V~	
Supply frequency range [Hz]	50/60 Hz	
Fuse designation	2 x T 3,15 AH, 250 V, UL-recognized	
Internal voltage supply	No	
Power consumption:	Current [A]	Power consumption [VA]
Normal operation 100 V/60 Hz	0.5 A	50 VA
Normal operation 240 V/50 Hz	0.3 A	72 VA
Max. allowed load of the additional socket / multiple sockets [A or VA]	1.6 A	
Protection class (I, II, III)	I	
Application part type (B, BF, CF):	BF (with Aquilex Fluid Control System inflow tube set) [1]	
Defibrillator protection (yes/no)	No	
Protection class (IP code)	IP41 (AQL-100P), IP21 (AQL-100CBS)	
Classification (I, IIa, IIb, III) in accordance with the Medical Device Directive 93/42/EEC, Appendix IX:	IIb	
Conformity with the following standards: (in the currently valid version)	IEC 60601-1 / EN 60601-1 / CAN/CSA-C22.2 NO. 60601-1 / AAMI ANSI ES 60601-1 IEC 60601-1-2 / EN 60601-1-2 / AAMI ANSI IEC 60601-1-2	
Operating conditions [°C] [°F], [%], [kPa]	10 to 40 °C/50 to 104 °F 30 to 70 % rel. humidity 70 to 106 kPa air pressure 3000 m max. altitude above sea level for device use	
Possible use with explosive anesthetic gases	This device is not designed for use with flammable anesthetic gases (Class AP) or flammable anesthetic gases with oxygen (Class APG).	
Storage conditions [°C] [°F], [%], [kPa]	5 to 40 °C/41 to 140 °F 5 to 85 % rel. humidity 70 to 106 kPa air pressure	
Transport conditions [°C] [°F], [%], [kPa]	-20 to 70 °C/-4 to 158 °F 5 to 90 % rel. humidity at 30 °C/86 °F 70 to 106 kPa air pressure	
Max. sound level [dB]	< 80 dB(A) (with acoustic signals)	
Max. load of container scale	25 kg/55,1 lbs, 4 containers	
Max. load of bag scale	12 kg/ 26.5 lbs, 2 bags 5 l each	
Maximum flow rate [l/min]	0.8 l/min ± 10 %	
Adjustable values		
	Pressure range [mmHg]	40-150 mmHg, 5.3-20 kPa
	Deficit limit [l]	0.6 - 2.5 l

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Measurement range		
	Pressure [mmHg]	0-500 mmHg
	Deficit [ml]	-995/+9995 ml
Accuracy		
	Pressure [mmHg]	± 10 mmHg
	Deficit [ml]	< 1 l: ± 60 ml > 1 l: ± 6 %
	Total inflow volume [ml]	± 10 %
Dimensions	Width x Height x Depth [mm], [in]	300 mm x 140 mm x 300 mm/11.8 in x 5.5 in x 11.8 in (AQL-100P), 670 mm x 1390 mm x 670 mm/26.4 in x 54.8 in x 26.4 in (AQL-100CBS)
Weight [kg], [lbs]		5.8 kg, 12.8 lbs (AQL-100P) 14.2 kg, 31.3 lbs (AQL-100CBS)
Interfaces:		
	ON/OFF signal for components	Irrigation Pump Unit: 1 x scale connection (flanged socket/round connection socket with 5 pins) 1 x service port (RS232 socket DSUB9/D-SUB9) Fluid Monitoring Unit/Container Scale: 1 x data connection (flanged female / round socket with 5 pins) for connection to bag scale. Fluid Monitoring Unit/Bag Scale: 2 x data ports (flanged socket/round connection socket with 5 pins): Connector to connect to pump, socket to connect with container scale
	Mains power socket	IEC 60320-1 C14
	RFID Transponder Technology	Working frequency: 13.5609 MHz, Transmitting power: -7.51 dBµA/m at 10 m
Essential Performance		Pressure build-up in the body cavity, control and measurement, limit: max. 150 mmHg (Normal Condition), 200 mmHg for max. 5 seconds (Single Fault Condition). Suction pressure: Limit value: min. - 84 kPa

[1] A tube set is not an application part in terms of the standard. However, it meets all the technical requirements for an application part.

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en

12 Guidelines and Manufacturer's Statement - Electromagnetic Compatibility

Precautionary measures

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

This device is to be used only for the purposes described in the manual and has to be installed, set up, and operated in compliance with the EMC notes and instructions.


WARNING!
The Aquilex Fluid Control System should not be used directly next to other devices as this could result in malfunctions. The Aquilex Fluid Control System was tested for compliance with IEC 60601-1-2 as a stand alone system. Therefore, do not stack other devices (e.g. MyoSure® Control Unit) on the system or the Irrigation Pump Unit. In particular, do not place any other device than the AQL-100P on the trays of the AQL-100CBS. If usage in the manner described above is nevertheless required, this system and the other devices should be monitored to make sure they function properly.

CAUTION!
Accessories
To ensure compliance with the requirements of IEC 60601-1-2 in the current version, the device Aquilex Fluid Control System must be used only with the accessories listed in chapter 13, Accessory List.

To ensure the basic safety and essential functionality in relation to electromagnetic interference over the life of the device, the device must be restarted after 24 hours so that a diagnostic self-test can be performed. The maintenance intervals indicated in Chapter 8.2 must also be observed.

12.1 Electrical Connections

Do not touch electrical connections identified with this warning label. Do not establish a connection between these plugs and sockets without first implementing precautionary ESD (electrostatic discharge) measures.

ESD (Electrostatic Discharge) precautionary measures

The following are ESD precautionary measures:

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

Hospital employees should be informed about and trained in ESD precautionary measures.

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12.2 Guidelines and Manufacturer's Statement – Electromagnetic Emissions

The Aquilex Fluid Control System is intended for use in the electromagnetic environment specified below. The user/operator of the Aquilex Fluid Control System should make sure the device is operated within such an environment.

WARNING!
Electromagnetic emissions may increase and rise above the permissible limits if other equipment (e.g. MyoSure® Control Unit) is stacked onto or placed directly next to the Aquilex Fluid Control System. The user is responsible for monitoring the devices to make sure they function properly.

Emitted interference measurements	Compliance	Electromagnetic environment guidelines
HF emission according to CISPR 11	Group 1	The Aquilex Fluid Control System uses HF energy solely for its internal functions. Therefore the HF emission is very low and it is unlikely that devices in close proximity will experience interference.
HF emission according to CISPR 11	Class B	The Aquilex Fluid Control System is suitable for use in all facilities including those in residential areas and those directly connected to a public utility network supplying buildings used for residential purposes as well.
Emission of harmonic oscillations according to IEC 61000-3-2	Class A	
Emission of voltage fluctuations / flickers according to IEC 61000 3-3	In compliance	

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12.3 Guidelines and Manufacturer's Statement/Electromagnetic Interference Immunity

TheAquilex Fluid Control System is intended for use in an electromagnetic environment as described below. The user/operator of the Aquilex Fluid Control System must make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	**IEC 60601-1-2 test level**	**Compliance Levels**	**Electromagnetic environment/guidelines**
Discharge of static electricity (ESD) according to IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be made from wood or concrete or covered with ceramic tiles. If the floor covering consists of synthetic material, the relative humidity should be at least 30%.
Electrical fast transients/bursts according to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines Modulation 100 KHz	± 2 kV for power lines ± 1 kV for input and output lines Modulation 100 KHz	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Voltage surges according to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Blackouts, brownouts, and fluctuations of the power supply according to IEC 61000-4-11	< 5% U_T* (> 95% dip in the U_T) for ½ cycle < 40% U_T (> 60% dip in the U_T) for 5 cycles < 70% U_T (> 30% dip in the U_T) for 25 cycles < 5% U_T (> 95% dip in the U_T) for 5 cycles	< 5% U_T * (> 95% dip in the UT) for 1/2 cycle 40% U_T (60% dip in the U_T) for 5 cycles 70% U_T (30% dip in the U_T) for 25 cycles < 5% U_T (> 95% dip in the U_T) for 5 seconds 0% U_T; 1/2 cycle at 0°, 45°, 90°, 135°, 180°, 225°,270°, and 315°0% U_T; 1 cycle and 70% U_T; 25/30 cycles Single phase: at 0°0% U_T; 250/300 cycles	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment. If the user/operator of the system requires the continuation of functionality after power interruptions/disruptions, it is recommended to supply the device with power from an uninterruptible power supply.
Supply frequency magnetic field (50/60 Hz) according to IEC 61000 4-8	30 A/m	30 A/m	Magnetic fields of the mains power frequency should comply with the typical values of business and hospital environments.

*Note: U_T is the mains alternating voltage before applying the test levels.

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12.4 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity

Electromagnetic interference immunity tests	Test level	Compliance	Electromagnetic environment guidelines
Conducted HF interference quantities according to IEC 61000-4-6 Radiated HF interference quantities according to IEC 61000-4-3	3 V_{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz	3 V_{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz	Portable and mobile wireless devices should not be used in closer proximity to the Aquilex Fluid Control System (including cables/lines) than the recommended safety distance calculated based on the transmitting frequency and the applicable formula. Recommended safety distance: $d = 1.2\sqrt{P}$ for 150 KHz to 80 MHz $d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ for 800 MHz to 2.7 GHz With **P** as the rated output of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and **d** as recommended safety distance in meters [m]. The field strength of stationary transmitters for all frequencies tested on site [a] should be lower than the concordance level. [b] Interference is possible in the proximity of devices featuring the following pictograph. 

Note 1: The higher frequency range applies for 80 and 800 MHz.

Note 2: These guidelines are probably not realizable in all cases. The distribution and spread of electromagnetic quantities differs depending on the absorption and reflection of buildings, objects, and people.

[a] The field strength of stationary transmitters such as base stations of wireless phones and cell phones, ham radio operators, AM and FM radio and TV stations can theoretically not always determined in advance. A study of the installation site should be considered to determine the electromagnetic environment concerning the stationary transmitter. If the measured field strength at the proposed Aquilex Fluid Control System installation and operation site exceeds the concordance levels listed above, the Aquilex Fluid Control System should be monitored to document proper functionality and operation as intended. If unusual performance characteristics are observed, additional measures may be required such as changing orientation or the location of the Aquilex Fluid Control System.

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[b] The field strength should be less than 3 V/m for the frequency range of 150 kHz to 80 MHz.

12.5 Recommended Safety Distances Between Portable and Mobile RF Telecommunications Devices and Aquilex Fluid Control System

The Aquilex Fluid Control System is intended for use in an electromagnetic environment where HF interferences are controlled. The user/operator of the Aquilex Fluid Control System can contribute to lowering electromagnetic emissions by complying with the minimum distance between portable and mobile HF telecommunications devices (transmitters) and the Aquilex Fluid Control System, depending on the output power of the communication equipment listed below.

Rated output of the transmitter [W]	Safety distance based on the transmitting frequency [m]		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	**80 MHz to 800 MHz** $d = 1.2\sqrt{P}$	**800 MHz to 2.7 GHz** $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The safety distance d in meters [m] for transmitters with a max. rated output not listed in the table above can be calculated by applying the corresponding formula in the respective column. P is the max. rated output of the transmitter in watts [W] as specified by the transmitter manufacturer.

Note 1: The higher frequency range applies at 80 MHz and 800 MHz.

Note 2: These guidelines are probably not realizable in all cases. The distribution and spread of electromagnetic quantities differs depending on the absorption and reflection of buildings, objects, and people.

WARNING!
Portable HF communication equipment can affect the performance characteristics of the Aquilex Fluid Control System. Such equipment must therefore comply with a minimum distance of 30 cm (regardless of all calculations) from the Aquilex Fluid Control System, its components and cables.

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13 Accessory List

The following accessories and peripherals are available:

Article	Order No.
Aquilex Fluid Control System complete tube set (inflow and outflow), disposable, sterilized using ethylene oxide	AQL-112
Aquilex Fluid Control System tube set vacuum (high and low)	AQL-114
Aquilex Fluid Control System container rings	AQL-200
Aquilex Fluid Control System power cord (US)	AQL-215
Aquilex Fluid Control System power cord (UK)	AQL-216
Aquilex Fluid Control System power cord (EU)	AQL-217

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
14 Warranty Information

Hologic warrants to the original purchaser of the Aquilex Fluid Control System that it shall be free of defects in material and workmanship when used as intended under normal surgical conditions and in conformance with its instructions for use and maintenance instructions. The obligation of Hologic under this warranty shall be limited to the repair or replacement, each at no charge, at the option of Hologic within one year from the date of purchase. Alternatively, Hologic may elect to repay or credit the original purchaser an amount equal to the purchase price of the defective equipment.

THIS WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES EXPRESSED OR IMPLIED INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR USE AND ALL OTHER OBLIGATIONS AND LIABILITIES ON THE PART OF HOLOGIC. HOLOGIC'S ENTIRE WARRANTY RESPONSIBILITY IS EXPRESSLY LIMITED TO REPAIR OR REPLACEMENT (AT HOLOGIC'S OPTION AND IN THE FORM ORIGINALLY SHIPPED) OF PRODUCT OR CORRECTION OF SERVICE SUBJECT TO ANY CLAIM, OR, AT HOLOGIC'S ELECTION, REPAYMENT OF, OR CREDITING CUSTOMER WITH, AN AMOUNT EQUAL TO THE HOLOGIC PRICE, FEE OR CHARGE THEREFOR. SUCH LIMITED WARRANTY IS GIVEN SOLELY TO THE ORIGINAL PURCHASER AND IS NOT GIVEN TO, NOR MAY IT BE RELIED UPON BY, ANY THIRD PARTY, INCLUDING, WITHOUT LIMITATION, CUSTOMERS OF PURCHASER. THIS WARRANTY IS VOID UPON TRANSFER OF PRODUCT BY PURCHASER TO ANY ENTITY WHO HAS LESS THAN FIFTY (50) PERCENT OWNERSHIP IN THE PRODUCT. THIS WARRANTY SHALL NOT APPLY TO AN AQUILEX SYSTEM OR TO THE AQUILEX FLUID CONTROL SYSTEM WHICH HAS BEEN SUBJECT TO ACCIDENT, NEGLIGENCE, ALTERATION, ABUSE, OR MISUSE, OR THAT HAS BEEN REPAIRED, MOVED, OR ALTERED BY ANYONE OTHER THAN AN AUTHORIZED HOLOGIC SERVICE PERSON. HOLOGIC MAKES NO WARRANTY WHATSOEVER WITH REGARD TO ACCESSORIES OR PARTS USED IN CONJUNCTION WITH THE AQUILEX FLUID CONTROL SYSTEM NOT SUPPLIED AND/OR MANUFACTURED BY HOLOGIC. THE TERM “ORIGINAL PURCHASER”, AS USED IN THE WARRANTY, SHALL BE DEEMED TO MEAN THAT PERSON OR ORGANIZATION AND ITS EMPLOYEES, IF APPLICABLE, TO WHOM THE AQUILEX SYSTEM WAS SOLD BY HOLOGIC.

Technical Support and Product Return Information

Contact Hologic Technical Support if the Aquilex Fluid Control System fails to operate as intended. If product is to be returned to Hologic for any reason, Technical Support will issue a Returned Materials Authorization (RMA) number. Return Aquilex System according to the instructions provided by Technical Support. Be sure to clean the Aquilex System with a clean damp cloth and germicide or isopropyl alcohol before returning it and include all accessories in the box with the returned unit.

Hologic and its distributors and customers in the European Community are required to comply with the Waste Electrical and Electronic Equipment (WEEE) Directive (2012/19/EC). Hologic is dedicated to meeting country specific requirements related to the environmentally sound treatment of its products. Hologic’s objective is to reduce the waste resulting from the disposal of its electrical and electronic equipment. Hologic realizes the benefits of subjecting such WEEE to potential reuse, treatment, recycling or recovery to minimize the amount of hazardous substances entering the environment. Hologic customers in the European Community are responsible for ensuring that medical devices marked with the following symbol, indicating that the WEEE Directive applies, are not placed into a municipal waste system unless authorized to do so by local authorities.

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Contact Hologic Technical Support to arrange for proper disposal of the Aquilex System in accordance with the WEEE Directive.

Hologic Technical Support

United States and Canada:

Phone: 1 800 442 9892 (toll-free) or 1 508 263 2900

Fax: 1 508 229 2795

Authorized European Representative:

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15 Glossary

Term	Statement
Embolism	Sudden capillary blockage due to embolus
Flow rate	Quantity (in ml) of irrigation fluid flowing through tube set per minute
Hypervolemia	An increased volume of circulating blood
Hyponatremia	A low concentration (< 130 mmol/l) of sodium in the patient's bloodstream
Hysteroscope	Endoscope to look inside the uterus
Intrauterine pressure	Pressure in uterine cavity
Intravasation	Entry of foreign matter into a blood vessel
Contamination Soiling	Pollution of rooms, water, foods, objects, or persons due to microorganisms or radioactive materials, biological poisons or chemical agents
Contraindication	Circumstances (e.g., age, pregnancy, certain illness, medication) prohibiting the use of an otherwise indicated measure (contrary to an indication)
Saline	Isotonic saline solution, i.e., one liter (l) contains 9.0 grams of sodium chloride.

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16 Appendix

16.1 Test Log

Date	Tests Performed	Results	Comment	Signature

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A
Authorized service technician 4

B
Bag change during surgery 28

C
Care and maintenance 4, 4
Certification 32
Connection to the wall outlet 13
Connect the fluid bags 26
Container change during surgery 28
Contamination 4
Contraindications 5

D
Deficit rate >300 ml/min 31

E
ESD (Electrostatic Discharge) precautionary measures 44
Exclusion of liability 4

G
Grounding contact 13

I
Initial system set-up 13
Inserting the tube set 26
Inspection tests 34
Instrument change during surgery 29
Intended Use 5

L
Loading/unloading scale while in operation 31

M
Manufacturer's specification 34

N
Normal Use 4

O
Only for U.S. operators 12, 13
Open outer packaging 25

P
Potential equalization 13
Precautionary measures 14, 14
Precise balancing 17
Pressure measuring and regulating 12
Pressure setting at restart 31

S
Safety threshold 26
Scale capacity 17
Scale overload 31
Serious system defect 31
Subject to technical changes 3

T
Technical application scope of the system 12
Testing the Flow Rate 35
Tube set recognition 21

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Deze gebruiksaanwijzing bevat auteursrechtelijk beschermde informatie waar copyright op bestaat. Alle rechten voorbehouden. Het is verboden om deze gebruiksaanwijzing zonder uitdrukkelijke, schriftelijke toestemming van W.O.M. WORLD OF MEDICINE GmbH geheel of gedeeltelijk door middel van fotokopieën, microfilm of met andere middelen te vermenigvuldigen of te verspreiden.